Smartline

RI Detector 2300/2400

Compact Value

The refractive index detectors **Smartline RI Detector 2300 and 2400** represent consistent advances in comparison to their very successful predecessors.

By refining the optical bench, crucial parameters such as drift, noise, and temperature stability were noticeably improved.

The design of both Smartline refractive index detectors allows for their use as stand-alone instruments as well as integrated into a Smartline Tower. The eluent connectors are easily accessible as they are located behind the front panel, which also includes the display and keypad. The redesigned control elements and the enlarged display markedly simplify the operation of these detectors.

The Smartline RI detectors can also be easily integrated into any manufacturer's system by utilization of the KNAUER open communication protocol. The connection is flexibly realized by means of the RS-232 interface or the remote control socket.

The KNAUER RI detectors of the Smartline series fulfill all the needs of modern refractive index detection systems. They are best suited for the detection of substances that do not absorb in the UV range, both in HPLC and gel permeation chromatography (GPC/SEC) applications.

The **Smartline RI Detector 2300** can be used in analytical applications for flow rates up to 5 ml/min. The **Smartline RI Detector 2400** has been developed for flow rates of up to 100 ml/min for semi-preparative use.

Digital data acquisition and control of the auto-zero feature, as well as rinsing functions, can be easily carried out from a PC. An additional analog output also makes it very simple to connect the detectors to an integrator or chart recorder.

Smartline RI Detector 2300/2400 – universal detection and an interesting alternative.

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
Technical Data and Ordering Information

Detection		
Detector Type	refractometer	
Light Source	LED	
Wavelength	950 nm	
Refractive Index Range	1.00–1.75 RIU	
Model	RI Detector 2300	RI Detector 2400
Sensitivity	3×10^{-8} RIU	4×10^{-7} RI
Noise	$\leq \pm 1.5 \times 10^{-8}$ RIU	$\leq \pm 2 \times 10^{-7}$ RIU
Measurement Range	$\pm 1 \times 10^{-3}$ RIU	$\pm 2 \times 10^{-3}$ RIU
Maximum Flow Rate	5 ml/min	100 ml/min
Measurement Angle	45°	15°
Cell Volume	15 μ l	9 μ l
Time Constants (both models)	0.1/0.2/0.5/1.0/2.0/5.0/10.0 s	

Technical Parameters

Analog Output (Integrator)	± 1 V (scalable in 16 increments)
Control	RS-232, remote control strip, touchpad
GLP Report	total operation time, serial number
Display LCD	2 lines, 24 characters each
Power Supply voltage range	90–260 V, 47–63 Hz
Dimensions	226 x 135 x 390 mm (W x H x D)
Weight	8.0 kg

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Order No.	Spare Parts and Accessories
A0287	15° Measurement Cell*
A0294	45° Measurement Cell*
A0277	Set of Seals*

* We recommend that exchange of cell or seals be done by authorized service personnel only.

Technical data are subject to change without notice.

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Visit www.knauer.net for details on complete HPLC systems, HPLC columns, and osmometers.

KNAUER · ASI · Advanced Scientific Instruments
Wissenschaftliche Gerätebau
Dr. Ing. Herbert Knauer GmbH
Hegauer Weg 38 · D-14163 Berlin, Germany
Telephone: +49 (0)30 / 80 97 27-0
Telefax: +49 (0)30 / 8 01 50 10
E-Mail: info@knauer.net
Internet: www.knauer.net

Your local distributor: