



KNAUER's monitoring and control software for AZURA® LC devices and systems

The hand-held Mobile Control (Chrom) allows a complete overview of all devices of the AZURA systems on one screen. Remotely check important parameters or control and monitor devices. The touch screen interface facilitates navigation using just your fingers. Choose Mobile Control as a basic, easy-to-use and cost-effective software solution!

Mobile Control provides full access to AZURA devices. Change device settings, set operating parameters, automate device control or check the system status and GLP data... Mobile Control features all functionalities of a device display.

Do you want to acquire data without the overhead of an advanced chromatographic data system?

Mobile Control Chrom features data acquisition from AZURA detectors in addition to full device control.

Help



System Overview in Mobile Control (Chrom)

Why Mobile Control software?

A software solution as device display, for device control and basic LC system operation.

Monitoring and Control So...



Help

Only pay for what you use

Mobile Control features basic functions to operate AZURA devices and systems. Mobile Control can operate dedicated applications which do not require a highly developed and cost-intensive Chromatographic Data System (CDS). Examples are GPC cleanup or other time-based purifications.

Save space

Mobile Control runs on a tablet. Especially in labs with little space avoiding a desktop PC with keyboard and monitor can be a decisive factor. The touch-optimized user interface

allows device control using just your fingers.

Save time

Mobile Control convinces due to an intuitive user interface and a clearly structured menu function. The training period is minimal in comparison to a complex CDS. Especially users with little experience in LC will easily get familiar with the software.

Free updates

With every release new features are available in Mobile Control. You can download the current version for free.

Free trial

To evaluate if Mobile Control holds up to your expectations, you can download Mobile Control and test the free trial option. Perfect for those who'd like to try before they buy.

Customized software design

Mobile Control is made by KNAUER and can be adapted to the requirements of our OEM partners.

[Download Mobile Control](#)

What is new in version 5.5.36?

We have released the new version of our award-winning software Mobile Control.

Convince yourself and test our full-featured trial version now.

[Help](#)

- P 6.1L isocratic with solvent selection valve supported
- Interactive queue and history for better accessibility of data
- Do not stop pump in case of disconnection of devices
- AS 6.1L - support of 100 µl syringe volume
- UVD 2.1S - support of further analog out scales
- New devices: Valve drive VU 4.1 and Assistant ASM 2.2L supported
- AS 6.1L support 30 vials tray
- Interactive queue and history for better accessibility of data
- Duplicate input facilitates the virtual keyboard entry
- Show gradient composition as varicolored background of chromatogram
- Optional second y-axis of chromatogram in Mobile Control & Data Viewer
- Select multiple programs for editing, exporting or deleting
- Support of Autosampler as standalone device/ Allow to start programs with first line only
- Enable autozero and valve switching during program run
- Allow instrument Power Up and stop pumps in Monitor Mode
- general improvements of user interface
- Sort program by date/name
- Redesigned reports
- Reports are saved in Results folder
- Allow flipped landscape rotation

Help





Control & Check
your system
remotely



Touch-optimized
user-interface




Option for data
acquisition



Free updates

Help

AVAILABLE PRODUCTS



Intuitive and cost-effective software for automated quality control of HPLC/FPLC columns
Mobile Control Chrom with data acquisition and column test option

Article No.: A9613



Display software for AZURA devices without data acquisition incl. tablet
Mobile Control without data acquisition including tablet and mount


Article No.: A9607



Display software for AZURA devices with data acquisition incl. tablet
Mobile Control Chrom with data acquisition including tablet and mount


Article No.: A9608

Help




**Display software for AZURA
devices without data acquisition**
Mobile Control without data
acquisition

Article No.: A9610



**Display software for AZURA
devices with data acquisition**
Mobile Control Chrom with data
acquisition

Article No.: A9612



**Upgrading Mobile Control to Mobile
Control Chrom**
Upgrading Mobile Control to Mobile
Control Chrom gaining data
acquisition

Article No.: A9614

Help



Tablet lock with stand
Tablet lock with stand for all tablets
Article No.: A9615



Tablet lock
Tablet lock for all tablets
Article No.: A9616



USB-LAN Adapter
USB-LAN ADAPTER Network
adapter USB 2.0 <-> 10/100
Ethernet for tablets
Article No.: A96181

DOWNLOADS

Help [ad installers of Mobile Control, Data Viewer & Firmware Wizard](#)

English

Instructions Mobile Control (Chrom) Version 5.5 (V6851-2)

Manual English

Release Note Mobile Control V5.5.35

Information English

Installation Information Mobile Control V5.5 (V6852A)

Information English

Update Installation Information Mobile Control V5.5 (V6853)

Information English

Update Installation Information Mobile Control V5.5 (V6853)

Information English

Upgrade Installation Information Mobile Control V5.5 (V6854)

Help ion English

Video tutorials

View all tutorials of KNAUER's award winning Mobile Control software.

[View tutorials](#)



We are happy to support you

Any questions about our products, solutions and services?

Need a quote?



+49 30 809727 0
Monday to Friday
7 am - 5 pm CET



sales@knauer.net
24 hours

Help